

Purification and treatment of recombinant human RyR2

MM Marco C. Miotto AM Andrew R. Marks

Updated date: Apr 27, 2023



An abbreviated version of this protocol was published in Science Advances in Jul 2022

Structural analyses of human ryanodine receptor type 2 channels reveal the mechanisms for sudden cardiac death and treatment

DOI: 10.1126/sciadv.abo1272

Detailed protocol

Homogenization buffer (H) 100 mL:

10 mM TRIS malate pH 6.8 (2 mL 0.5 M), 1 mM EGTA (200 μ L 0.5M), DTT 1 mM (100 μ L 1 M or 15.4 mg), 1 mM Benzamidine (18 mg), 0.5 mM AEBSF (12 mg), 2 protease inhibitor (PI) tablet. Add GST-Cal2.

Buffer S 15 mL:

10 mM HEPES pH 7.4 (300 μ L 0.5 M), 0.85 M NaCl (0.75 g), 1.5% CHAPS (225 mg), 0.5% PC (7 μ L, stock: 10% in 10% CHAPS), 1 mM EGTA (3 μ L 0.5M), 2 mM DTT (30 μ L 1 M or 4.63 mg), 0.5 mM AEBSF (6 mg), 1 mM Benzamidine (8 mg), 1 protease inhibitor (PI) mini tablet. Add GST-Cal2.

Buffer D 50 mL:

10 mM HEPES pH 7.4 (1000 μ L 0.5 M), 1.0% CHAPS (500 mg), 0.1% PC (0.5 mL, stock: 10% in 10% CHAPS), 1 mM EGTA (100 μ L 0.5M), 2 mM DTT (1 μ L 1 M), 1 protease tablet.

Buffer Wa 500 mL:

10 mM HEPES pH 7.4 (10 mL 0.5 M), 0.4% CHAPS (2 g), 1 mM EGTA (1 mL 0.5M), 0.5 mM TCEP (65 mg), 0.001% DOPC (evaporate 0.2 mL from chloroform solvent in glass tube and dissolve in buffer W), 230 mM NaCl (6.7 g).

Buffer Wb (100 mL Wa):

NaCl 666 mM (+2.5 g).

Buffer Ge (100mL Wa):

GSH 10mM (0.3g), DTT 1mM (0.1mL 1M), check pH 7.7-8.0.

Protocol:

1. Resuspend cell pellet in buffer H, approx: 1.5 mL/dish. Add GST-Cal2 (40 nmoles, 100 equivalents of 1mg of RyR2). Incubate 15 min. Lysate via sonication (50%, pulses: 2 min; 15 s on, 20 s off).
2. Centrifuge to separate DNA and debris: 10' 3000 xg.
3. Centrifuge to separate membranes: 30' 100.000 xg (30k rpm Ti-45, 34k rpm Ti-55 or Ti-70).
4. Resuspend pellet in 10-15 mL S, adding GST-Cal2 (40 nmoles). Homogenize.
5. Add 35-50 mL D with glass homogenizer. Incubate for 30-60 min in overhead shaker (cold room).
6. Centrifuge to separate insoluble debris: 30' 100.000 xg.
7. Vacuum filter and load onto Q column (5mL), pre-equilibrated with Wa, at 1 mL/min
8. Wash with Wa for 5 CV
9. Elute 3 volumes with isocratic buffer at 40 mS/cm (60-70% buffer Wb) and pool RyR2 fractions.
10. Add 80 nmoles of GST-Cal2.
11. Load into GSTrap column, pre-equilibrated with Wa, at 0.2-0.5 mL/min. Leave ON recirculating.
12. Wash with Wa₂₀₀ for 5-10 CV, at 1-2 mL/min. Attach Q column (1mL). Wash with Wa₂₀₀ for 5 CV
13. Elute with Ge for 5-10 CV or until stabilization of UV. Wash with Wa.
14. Elute with a linear gradient between 30-45 mS/cm NaCl (40-70% buffer Wb). Pool RyR2 fractions.
15. Thrombin digestion: add 50 U of Thrombin. PKA phosphorylation: 100 U PKA, 10 mM EGTA (60 μ L 0.5 M), 8 mM MgCl₂ (24 μ L 1 M), 100 μ M ATP (1.6 μ L 200 mM). Dephosphorylation: 6-12 μ L phosphatase lambda, buffer 10X, Mn 10X. Incubate 30 min at RT.
16. Concentrate and load into SEC (TSKgel). Pool RyR2 fractions.
17. Concentrate to 20 μ L and filter (centrifugal filters 0.22 μ m).
18. Measure concentration by NanoDrop.

How to cite:(Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Miotto, M. and Marks, A. (2023). Purification and treatment of recombinant human RyR2. Bio-protocol Preprint. bio-protocol.org/prep2254.
2. Miotto, M. C., Weninger, G., Dridi, H., Yuan, Q., Liu, Y., Wronska, A., Melville, Z., Sittenfeld, L., Reiken, S. and Marks, A. R.(2022). Structural analyses of human ryanodine receptor type 2 channels reveal the mechanisms for sudden cardiac death and treatment. Science Advances 8(29). DOI: [10.1126/sciadv.abo1272](https://doi.org/10.1126/sciadv.abo1272)

Copyright: © 2023 The Author(s); This is an open-access article under the CC BY license (<https://creativecommons.org/licenses/by/4.0/>).